


Disposable Laparoscopic Hand/Foot Control Instrument

Medstar's Laparoscopic electrosurgical instruments provide excellent dissection and coagulation for a variety of laparoscopic procedures. The electrodes are made from 33 cm and 45 cm solid stainless steel shaft. Safe delivery of electrosurgery with PTFE and Polyolefin double insulation reduces the potential of unintended site burns. 45 cm length for bariatric procedures, various tip configurations offer options for range of laparoscopic procedures. Compatible with 5 mm standard laparoscopic ports. Compatible with standard hand switching pencils, also compatible with standard foot control cord, electrode interchangeable in the sterile field.



Disposable Laparoscopic Electrode

Laparoscopic Electrode Compatible with Foot Control Cord



Laparoscopic Electrode With Banana Plug Compatible with Foot Control Cord		330mm Length	450mm Length
	L-Wire Teflon Coated	EES01-LWR330G	EES01-LWR450G
	L-Hook Teflon Coated	EES01-LHK330G	EES01-LHK450G
	J-Wire Teflon Coated	EES01-JWR330G	EES01-JWR450G
	J-Hook Teflon Coated	EES01-JHK330G	EES01-JHK450G
	O-Wire Teflon Coated	EES01-OWR330G	EES01-OWR450G
	Straight Blade Teflon Coated	EES01-SBL330G	EES01-SBL450G
	Curved Blade Teflon Coated	EES01-CBL330G	EES01-CBL450G
	Spatula Teflon Coated	EES01-SP330G	EES01-SP450G
	Curved Spatula Teflon Coated	EES01-CSP330G	EES01-CSP450G
	Needle Teflon Coated	EES01-NDL330G	EES01-NDL450G
	Ball Teflon Coated	EES01-BL330G	EES01-BL450G

Laparoscopic Electrode Compatible with Standard Electrosurgical Pencil




Laparoscopic Electrode With Hexagonal Plug Compatible with Standard Electrosurgical Pencil		330mm Length	450mm Length
	L-Wire Teflon Coated	EES02-LWR330G	EES02-LWR450G
	L-Hook Teflon Coated	EES02-LHK330G	EES02-LHK450G
	J-Wire Teflon Coated	EES02-JWR330G	EES02-JWR450G
	J-Hook Teflon Coated	EES02-JHK330G	EES02-JHK450G
	O-Wire Teflon Coated	EES02-OWR330G	EES02-OWR450G
	Straight Blade Teflon Coated	EES02-SBL330G	EES02-SBL450G
	Curved Blade Teflon Coated	EES02-CBL330G	EES02-CBL450G
	Spatula Teflon Coated	EES02-SP330G	EES02-SP450G
	Curved Spatula Teflon Coated	EES02-CSP330G	EES02-CSP450G
	Needle Teflon Coated	EES02-NDL330G	EES02-NDL450G
	Ball Teflon Coated	EES02-BL330G	EES02-BL450G

Disposable Laparoscopic Hand/Foot Control Instrument

Laparoscopic Electrode with Foot Control Cord attached



Disposable Laparoscopic Electrode with Foot Control Cord attached (7.9mm pin plug)		330mm Length	450mm Length
	L-Wire Teflon Coated	EES01-LWR330GC1	EES01-LWR450GC1
	L-Hook Teflon Coated	EES01-LHK330GC1	EES01-LHK450GC1
	J-Wire Teflon Coated	EES01-JWR330GC1	EES01-JWR450GC1
	J-Hook Teflon Coated	EES01-JHK330GC1	EES01-JHK450GC1
	O-Wire Teflon Coated	EES01-OWR330GC1	EES01-OWR450GC1
	Straight Blade Teflon Coated	EES01-SBL330GC1	EES01-SBL450GC1
	Curved Blade Teflon Coated	EES01-CBL330GC1	EES01-CBL450GC1
	Spatula Teflon Coated	EES01-SP330GC1	EES01-SP450GC1
	Curved Spatula Teflon Coated	EES01-CSP330GC1	EES01-CSP450GC1
	Needle Teflon Coated	EES01-NDL330GC1	EES01-NDL450GC1
	Ball Teflon Coated	EES01-BL330GC1	EES01-BL450GC1



Reusable Laparoscopic Foot Control Cord
7.9mm pin plug

Hand Control Laparoscopic Electrodes



Hand Control Laparoscopic Electrodes (3-pin plug)		330mm Length	450mm Length
	L-Wire Teflon Coated	PL01-LWR330G	PL01-LWR450G
	L-Hook Teflon Coated	PL01-LHK330G	PL01-LHK450G
	J-Wire Teflon Coated	PL01-JWR330G	PL01-JWR450G
	J-Hook Teflon Coated	PL01-JHK330G	PL01-JHK450G
	O-Wire Teflon Coated	PL01-OWR330G	PL01-OWR450G
	Straight Blade Teflon Coated	PL01-SBL330G	PL01-SBL450G
	Curved Blade Teflon Coated	PL01-CBL330G	PL01-CBL450G
	Spatula Teflon Coated	PL01-SP330G	PL01-SP450G
	Curved Spatula Teflon Coated	PL01-CSP330G	PL01-CSP450G
	Needle Teflon Coated	PL01-NDL330G	PL01-NDL450G
	Ball Teflon Coated	PL01-BL330G	PL01-BL450G